



Brandon Medical®

Member of REINSBERG GROUP

Anteris™ Ultra Accessories



Bringing innovation to life



SMART Solutions for Healthcare Spaces

World-leading designer and manufacturer of capital medical equipment: operating theatres, critical care areas & primary care.

From standalone products to integrated solutions, our connected equipment portfolio empowers construction, specifier, architect and estates and facilities management teams.

Brandon Medical is a UK company that delivers SMART turnkey equipment for acute and primary healthcare worldwide. For 75 years, British engineering skills and a strong work ethic have fuelled continuous product innovation and development to provide healthcare professionals with reliable, high-quality, and affordable medical equipment packages for operating theatres and critical care in over 70 countries around the globe.

We are acknowledged experts in manufacturing medical lighting and control systems and designing medical supply pendants, medical power, and medical audio-video systems according to customers' actual needs.

We specialise in enhancing healthcare environments by combining knowledge-driven intuition, decades of expertise, and practical innovation supported by continuous research and development.

Our 360-degree portfolio supports clinical teams and informs future intelligence, enabling the implementation of advanced best practices for better patient outcomes.

Constantly evolving, we help design the hospitals of tomorrow and drive healthcare excellence.

Contents

- 04. General Accessories
- 14. Head Positioning
- 18. Arm Positioning
- 22. Leg Positioning
- 28. Neurosurgery
- 30. Orthopaedics and Traumatology
- 38. Positioning Aids
- 42. Accessories for Anteris Ultra 700 Series
- 50. Accessories for Anteris Ultra 600 Series

General Accessories



Rotatable clamp

Art.-no. 61511

- For fixation of accessories with round rod (18mm)
- Material: Stainless steel



Rotatable clamp with quick lock and release lever

Art.-no. 61512

- For fixation of accessories with round rod (18mm)
- Material: Stainless steel, Aluminium



Clamp for square rod (25 x 10mm)

Art.-no. 61551

- Material: Aluminium



Clamp for square rod (20 x 20mm)

Art.-no. 61552

- Material: Aluminium



Clamp for round rod (18mm)

Art.-no. 61553

- Material: Aluminium



**Clamp with quick lock and release lever
for round rod (18mm)**

Art.-no. 61554

- Material: Stainless steel, Aluminium

General Accessories



Anaesthesia screen

Art.-no. 61750

- Extendable to both sides
- Incl. clamp
- Material: Stainless steel



Arm restraint cuff

Art.-no. 61755

- For attachment at anaesthesia screen
- Suitable for: 61750
- Material: Stainless steel



Flexible anaesthesia screen

Art.-no. 61751

- Bendable
- Individually adjustable
- Incl. clamp
- L: 2,000mm
- Material: Stainless steel



Tube support

Art.-no. 61752

- For fixation of breathing tube
- Bendable support is infinitely variable
- Adjustable in height
- Incl. clamp
- Material: Stainless steel



Infusion holder

Art.-no. 61800

- Extendable
- Four hooks
- Incl. clamp
- Material: Stainless steel



Rotatable infusion holder

Art.-no. 61805

- Extendable
- Rotatable
- Four hooks
- Incl. rotatable clamp
- Material: Stainless steel

General Accessories



Mayo table

Art.-no. 61410

- Incl. clamp
- Adjustable in height, rotatable and inclinable
- L x W: 450 x 350mm
- Material: Stainless steel



Side rail extension

Art.-no. 61560

- Incl. clamp
- Extends side rail up to 320mm
- For attachment of accessories at side rail (EU)
- Material: Stainless steel



Sliding handles (pair)

Art.-no. 61571

- For easy and comfortable movement of operating table
- Incl. clamp
- Material: Stainless steel



Side support (pair)

Art.-no. 61341

- Incl. clamp
- Upholstery: PUR, L x W each 220 x 110mm, electrically conductive
- Material: Stainless steel



Shoulder support (pair)

Art.-no. 61351

- Incl. clamp
- Adjustable in height and extendable
- Upholstery: PUR, L x W each 220 x 110mm, electrically conductive
- Material: Stainless steel



Shoulder support incl. concave upholstery (pair)

Art.-no. 61352

- Incl. clamp
- Adjustable in height and extendable
- Upholstery: PUR, L x W each 220 x 110mm, electrically conductive
- Material: Stainless steel

General Accessories



Shoulder support incl. straight upholstery (pair)
Art.-no. 61353

- Incl. clamp
- Adjustable in height and extendable
- Upholstery: PUR, L x W each 220 x 110mm, electrically conductive
- Material: Stainless steel



Double side support
Art.-no. 61362

- Incl. clamp
- Support of pelvis section
- Adjustable in height and extendable
- Upholstery: PUR, Ø 100mm, electrically conductive
- Material: Stainless steel



Joint arm for body supports
Art.-no. 61367

- Flexible positioning of body supports by means of three joints
- Quick lock and release of joint arm by means of central screw
- Material: Stainless steel



Holder for body supports

Art.-no. 61361

- Adjustable in height, extendable and swivelling (360°)
- Positioning of body supports by means of ball joint
- Material: Stainless steel



Body support for pubis, sacrum and sternum

Art.-no. 61361_1

- Upholstery: PUR, Ø 100mm, electrically conductive
- Suitable for: 61361, 61367 and 81367
- Material: Stainless steel



Body support for side and shoulder

Art.-no. 61361_2

- To support side and shoulder
- Upholstery: PUR, L x W 220 x 110mm, electrically conductive
- Suitable for: 61361, 61367 and 81367
- Material: Stainless steel

General Accessories



Body support for back and buttocks

Art.-no. 61361_3

- Upholstery: PUR, L x W 220 x 110mm, electrically conductive
- Suitable for: 61361, 61367 and 81367
- Material: Stainless steel



Supporting roll

Art.-no. 61361_4

- For lateral positioning of patient
- Upholstery: PUR, L x Ø 190 x 100mm, electrically conductive
- Suitable for: 61361, 61367 and 81367
- Material: Stainless steel



Instrument disposal round form

Art.-no. 91611


- Incl. clamp
- Swivelling
- Ø 280mm, Material: Stainless steel



Instrument disposal squared form

Art.-no. 91612

- Incl. clamp
- Swivelling
- L x W x H: 325 x 265 x 65mm, Material: Stainless steel



Bowl with drainage and tube

Art.-no. 91613

- Incl. clamp
- Swivelling
- Ø 320mm
- Material (Bowl): Stainless steel



Bucket

Art.-no. 91608

- With lid and tube-inlet
- Volume: 12 litre
- Material: Stainless steel
- Suitable for: 91613



Hand switch support

Art.-no. 61575

- To place hand switch at lateral or head-sided side rail
- Tiltable
- Material: Stainless steel

Head Positioning



Universal adapter for Anteris Ultra 600 Series
Art.-no. 61100

- For attachment of a variety of head positioning and neurosurgical accessories
- Material: Stainless steel



Universal adapter for Anteris Ultra 700 Series
Art.-no. 71100

- For attachment of a variety of head positioning and neurosurgical accessories
- Material: Stainless steel



Horseshoe shaped head rest with cut-out
Art.-no. 61114_1

- For positioning head in prone or supine position
- Oval form
- Tiltable and rotatable by means of ball joint
- Adjustable in height
- Upholstery: PUR, electrically conductive
- Material: Stainless steel
- Attachment at 61100, 61385, 61390, 71100 or 71390



Horseshoe shaped head rest two-parted
Art.-no. 61114_3

- For positioning head in prone or supine position
- Incl. detachable upholstery and fixture strap
- Extendable in width
- Rotatable and tiltable by means of ball joint
- Upholstery: Viscoelastic foam, electrically conductive
- Material: Stainless steel
- Attachment at 61100, 61385, 61390, 71100 or 71390



Head support

Art.-no. 61395

- For comfortable and safe positioning of patient's head during shoulder surgeries
- Attachment at tripartite back rest or universal adapter
- Extendable
- Incl. upholstery



Joint arm

Art.-no. 61117

- For attachment of head positioning accessories
- Flexible positioning of body supports by means of three joints
- Quick lock and release of joint arm by means of central screw
- Attachment at universal adapter 61100 or 71100
- Material: Stainless steel



Head support for joint arm

Art.-no. 61117_2

- For stable positioning and fixation of patient's head
- Incl. upholstery and fixation strap
- Upholstery: Viscoelastic foam, electrically conductive
- Attachment at 61117

Head Positioning



Horseshoe shaped head rest with cut-out for joint arm
Art.-no. 61117_3

- For positioning head in prone or supine position
- Oval form
- With cut-out
- Upholstery: PUR, electrically conductive
- Attachment at 61117



Horseshoe shaped head rest two-parted for joint arm

Art.-no. 61117_5

- Incl. detachable upholstery and fixture strap
- For positioning head in prone or supine position
- Extendable in width
- Upholstery: Viscoelastic foam, electrically conductive
- Attachment at 61117



Wrist support ring for head positioning accessories

Art.-no. 61113

- With connecting piece for horseshoe shaped head rest
- Material: Stainless steel
- Suitable for: 61114_1, 61114_3, 61117_3, 61117_4, 61117_5



Arm Positioning



Arm rest

Art.-no. 61211

- Incl. clamp and fixture straps
- Rotatable and tiltable by means of ball joint
- Adjustable in height and horizontally slidable
- Upholstery: Viscoelastic foam, L x W x H: 400 x 115 x 80mm, electrically conductive, detachable
- Material: Stainless steel



Beach-Chair- shoulder OP arm positioner

Art.-no. 61227

- Incl. clamp and fixture straps
- For positioning of patient's arm in beach chair position
- Quick lock and release of joint arm by means of central screw
- Upholstery: Viscoelastic foam, L x W x H: 500 x 115 x 80mm, electrically conductive, detachable
- Material: Stainless steel



Arm rest with joint arm and one-hand control

Art.-no. 61225

- Incl. clamp and fixture straps
- Three joints enable flexible positioning of patient's arm
- Quick lock and release of joint arm by means of central screw
- Upholstery: Viscoelastic foam, L x W x H: 500 x 115 x 80mm, electrically conductive, detachable
- Material: Stainless steel



**Arm rest with one-hand control and comfort upholstery
(Length 600mm)**

Art.-no. 61222

- Incl. clamp and fixture straps
- Rotatable by 360°, swivelling and tiltable by means of one-hand control
- Adjustable in height
- Upholstery: Viscoelastic foam, L x W x H: 600 x 165 x 40mm, electrically conductive, detachable
- Material: Stainless steel



**Arm rest with one-hand control and comfort upholstery
(Length 450mm)**

Art.-no. 61223

- Incl. clamp and fixture straps
- Rotatable by 360°, swivelling and tiltable by means of one-hand control
- Adjustable in height
- Upholstery: Viscoelastic foam, L x W x H: 450 x 165 x 40mm, electrically conductive, detachable
- Material: Stainless steel



Arm rest standard

Art.-no. 61220SF

- Incl. clamp and fixture strap
- Swivelling
- Material: Stainless steel
- Upholstery: Viscoelastic foam, L x W x H: 600 x 165 x 60mm, electrically conductive
- Also available with PUR upholstery, H = 60mm (Art.-no. 61220)

Arm Positioning



Side support / arm rest

Art.-no. 61092

- For stabilisation of patient or positioning of arm
- Hinged in three steps (upwards, horizontally, downwards)
- L x H: 600 x 140 / 235mm
- Upholstery made of artificial leather covered foam, H= 20mm has to be ordered separately (art.-no. 56031)



Arm protection shield made of stainless steel

Art.-no. 71345

- For safe positioning and fixation of the patient's arm
- Easy cleaning due to removable upholstery
- Dimensions (W x H): 420 x 160mm



Hand cuffs (pair)

Art.-no. 61400

- Incl. clamp
- For fixation of hands or arms
- L x W: each 180 x 80mm
- Material: Stainless steel



Leg Positioning



Leg rest according to Goepel (piece)

Art.-no. 61700

- For positioning patient in lithotomy position
- Upholstery: PUR, electrically conductive
- Material: Stainless steel
- Rotatable clamp (61511 or 61512) and fixture straps (B81600M1) to be ordered separately



Leg rest according to Goepel Comfort (pair)

Art.-no. 81705

- For positioning patient in lithotomy position
- Upholstery: PUR, electrically conductive
- Material: Stainless steel
- Rotatable clamp (61511 or 61512) and fixture straps (B81600M1) to be ordered separately



Stirrups with one-hand-adjustment, gas spring assisted (pair)

Art.-no. 61725

- Incl. pad and silicone fastening straps (no Velcro straps)
- For positioning the patient in lithotomy position
- Stepless positioning of the legs
- Gas spring assisted one-hand operation
- Max. Patient weight: 159kg (350 lbs)
-

Art.-no. 61726

- Max. Patient weight: 227kg (500 lbs)

Art.-no. 61729

- Max. Patient weight: 363kg (800 lbs)

*2x art. no. 61727 or 2x art. no. 61728
are required for attachment to the OR table*



Secure-Release™ Clamp, clip-on, EU

Art.-no. 61727

- suitable for item-No. 61725, 61726, 61729



Select-Blade Clamp, suspendable, EU

Art.-no. 61728

- suitable for item-No. 61725, 61726, 61729



Leg rest for adiposity patients (pair)

Art.-no. 61710

- Incl. upholstery and fixture strap
- Especially for the positioning of adipose patients in lithotomy position
- Tiltable by means of ball joint
- Upholstery: 300 x 285 x 155 / 100mm, electrically conductive
- Material: Stainless steel
- Rotatable Clamp (61511 or 61512) to be ordered separately

Leg Positioning



Elbow support (pair)

Art.-no. 61780

- Incl. clamp
- To support elbows of surgeon during procedures in lithotomy position
- L x W: 145 x 120mm
- Upholstery: PUR, electrically conductive
- Material: Stainless steel



Lower leg fixation

Art.-no. 61760

- Incl. clamp
- Upholstery: PUR, electrically conductive
- Material: Stainless steel



Thigh fixation

Art.-no. 61770

- Incl. clamp
- Upholstery: PUR, electrically conductive
- Material: Stainless steel



**Rectal positioning device for Anteris Ultra 600 Series
Art.-no. 61632**

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Adjustable in height (200mm)
- Integrated side rails (EU)
- Safe working load: 150kg
- Upholstery (lower): L x W x H 370 x 500 x 490 - 690mm, electrically conductive
- Material: Stainless steel



**Rectal positioning device for Anteris Ultra 700 Series
Art.-no. 71632**

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Adjustable in height (200mm)
- Integrated side rails (EU)
- Safe working load: 150kg
- Upholstery (lower): L x W x H 370 x 500 x 490 - 690mm, electrically conductive
- Material: Stainless steel



**Knee supports for Anteris Ultra 600 Series (piece)
Art.-no. 61150**

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Adjustable in height and swivelling
- Upholstery: PUR, electrically conductive
- Material: Stainless steel

Leg Positioning



Rectoscopy roll

Art.-no. 61315

- Incl. clamp
- Adjustable in height
- Upholstery: PUR, electrically conductive
- Material: Stainless steel



Knee support

Art.-no. 61630

- For positioning patient in knee-elbow position due to proctologic treatments or spinal surgery
- Incl. clamp
- Upholstery: PUR, L x W 525 x 245mm, electrically conductive
- Material: Stainless steel



Foot plates (pair)

Art.-no. 61631

- Incl. clamps
- To stabilize patient in Anti-Trendelenburg position
- For extension of lying surface (foot plates can also be hinged down)
- Upholstery: Viscoelastic foam, L x W x H each 270 x 270 x 60mm, electrically conductive, detachable
- Material: Stainless steel



Body and leg strap

Art.-no. 61405

- Incl. clamp
- Steplessly adjustable
- L x W: 2,300 x 80mm
- Material: Stainless steel



Leg strap (pair)

Art.-no. 61406

- For fixation of legs at divided leg rest
- Steplessly adjustable



Leg strap (pair)

Art.-no. 61407

- Incl. clamp
- For divided leg rests
- L: each 1,200mm
- Safety load: 227kg

Neurosurgery



Standard holder

Art.-no. 661065

- For fixation of skull clamp or horseshoe headrest
- Adjustable with one hand
- Large adjustment radius
- Very low weight
- Only in connection with 61100 or 71100



Swivel adapter

Art.-no. 661067

- Connects standard holder (661065) with skull clamp (661069)
- Enables rotation adjustment of skull clamp
- Suitable for: 661065, 661069, 661076 and 661078



Skull clamp

Art.-no. 661069

- For fixation of head during neurosurgical procedures
- Locking mechanism enables quick and secure fixation
- Only in connection with 661065 and 661067
- Skull clamps (661073 and 661074) are not included in delivery



Skull pins adults

Art.-no. 661073

- Can be sterilised, for multiple uses, three pieces, stainless steel
- Suitable for 661069

Skull pins children

Art.-no. 661074

- Same as 661073 but for children



Cross bar adapter

Art.-no. 661071

- Incl. clamps
- For secure fixation of head in seated position of patient
- Adaptable to the width of operating table
- Suitable for: 661065



Horseshoe headrest for adults

Art.-no. 661076

- Incl. detachable gel pads
- For positioning of head in prone or supine position
- Only in connection with 661065 and 661067

Horseshoe headrest for children

Art.-no. 661078

- Same as 661076, but for children

Orthopaedics and Traumatology



Extension device for Anteris Ultra 600 Series Art.-no. 63000_1

- Levitating
- Max. patient weight: levitating 136kg, with supporting bars 250kg
- Material: Stainless steel
- Alternatively available with viscoelastic foam upholstery, H = 80mm (Art.-no. 63000_1SF)

Included in delivery:

- Fixation adapter for operating table incl. pelvis part with upholstery
- Countertraction post for perineum and for neck of femur
- Two swivelling extension bars with double joints
- Two swivelling and height adjustable drawing spindles with adaption for extension shoe holders
- One trolley with storage basket for additional accessories
- Two height adjustable supporting bars



Extension device for Anteris Ultra 700 Series Art.-no. 63000_2

- Same as 63000_1, but for Anteris Ultra 700



Coupling point for extension device 63000_2 for Anteris Ultra 700 Series operating tables without sliding
Art.-no. 63005_1

- For attachment of extension device to operating table

Coupling point for extension device 63000_2 for Anteris Ultra 700 Series operating tables with sliding
Art.-no. 63005_2

- For attachment of extension device to operating table



Extension shoe for adults (piece)
Art.-no. 61010

- For fixation of foot
- Suitable for: 63000_1, 63000_2

Extension shoe for children
Art.-no. 61020

- Same as 61010 but for children



Countertraction post for femur for extension device 63000_1
Art.-no. 63025_1

- For extension in lateral position
- Radiolucent perinneeal arc

Countertraction post for femur for extension device 63000_2
Art.-no. 63025_2

- Same as 63025_1 but for extension device 63000_2

Orthopaedics and Traumatology



Countertraction post for tibia and fibula
Art.-no. 63045

- For positioning and fixation of patient during surgery of lower leg
- Radiolucent padded roll
- Adjustable in height and inclinable
- Material: Stainless steel
- Suitable for: 63000_1, 63000_2



Pelvis plaster cast unit
Art.-no. 63030

- Usable for adults and children
- Slidable countertraction post
- Material: Stainless steel
- Suitable for: 63000_1, 63000_2



Horseshoe holder "Kirschner"
Art.-no. 63060

- For extensions using Steinmann nails
- Material: Stainless steel
- Suitable for: 63000_1, 63000_2



Leg support
Art.-no. 63026

- Incl. fixture straps
- For positioning of healthy leg during extension
- Tiltable upholstery by means of ball joint
- Upholstery: Viscoelastic foam, L x W x H: 400 x 115 x 80mm, electrically conductive, detachable
- Material: Stainless steel
- Suitable for: 63000_1, 63000_2



Transfer leg rest (pair)
Art.-no. 63070_1

- For preoperative positioning of legs
- Incl. upholstery, L x W x H: 600 x 200 x 35mm, electrically conductive
- Suitable for: 63000_1

Transfer leg rest
Art.-no. 63070_2

- Same as 63070_1, but for 63000_2



Trolley for extension device
Art.-no. 63090

- For storage and transport of extension device and accessories
- Material: Stainless steel
- Suitable for: 63000_1, 63000_2

Orthopaedics and Traumatology



Radiolucent extension device for Anteris Ultra 600 Series

Art.-no. 64000

Consisting of:

- Connecting plate made of carbon fibre with outer segments which can be detached and hinged down
- Countertraction post made of carbon fibre with upholstery roll \varnothing : 90mm / H: 220mm
- Swivelling leg rests incl. attachment for traction device with segments which can be detached and hinged down
- Traction device with rotary clamp for attachment of extension shoe
- Extension shoe (1 pc.) with mounting clamp and disposable inlay set (10 pcs.)
- Trolley for storage of extension device and accessories



Extension shoe

Art.-no. 64010

- (Additional) extension shoe for extension devices
- With mounting clamp
- For use with disposable inlay set 10 pcs. (Art.-no. 64011)



Arm extension device

Art.-no. 62050

- For intramedullary nailing of humerus
- Radiolucent countertraction post
- Incl. tension rod, countertraction post and clamp



Weinberger hand cuff

Art.-no. 62055

- Secure fixation of hand and upper arm during intramedullary nailing of humerus
- Radiolucent
- Can be sterilised



Upper arm support

Art.-no. 61226

- For positioning of upper arm
- Radiolucent
- Adjustable in height

Orthopaedics and Traumatology



Arm and hand surgical table levitating

Art.-no. 61231

- Incl. clamps and upholstery
- Adjustable in height and radiolucent
- Upholstery: Artificial leather covered foam, L x W x H: 600 x 340 x 35mm, electrically conductive
- Material: Stainless steel, plastic



Arm and hand surgical table with supporting leg

Art.-no. 61232

- Incl. clamp, upholstery and supporting leg
- Height adjustment from 700 to 1,150mm
- Table top made of radiolucent plastic, dimensions (L x W x H): 815 x 520 x 5mm
- Upholstery: SF, H = 30mm, electrically conductive
- Material: Stainless steel, plastic



Arm and hand surgical table made of carbon

Art.-no. 81230

- Incl. clamp, upholstery and supporting leg
- Height adjustment from 720 to 1,150mm
- Table top made of radiolucent carbon, dimensions (L x W x H): 820 x 525 x 30mm
- Upholstery: SF, H = 30mm, electrically conductive
- Material: Stainless steel



Meniscotomy roll

Art.-no. 61330

- Incl. clamp
- Adjustable in height
- Upholstery: PUR, electrically conductive
- Material: Stainless steel



Arthroscopic leg holder system

Art.-no. 61328

- Incl. upholstery, clamp and fixture strap
- For fixation of thigh
- Horizontally and vertically adjustable
- Width of leg holder adjustable in three steps
- Upholstery: Artificial leather covered foam, electrically conductive
- Material: Stainless steel



Knee positioning device

Art.-no. 61327

- Incl. upholstery, clamps and fixture strap
- Intraoperative flexion of knee and fixation of thigh
- Electrically adjustable by means of foot switch
- Width of leg holder adjustable in three steps
- Suitable for different leg lengths
- Tilt angle: 0° - 60°
- Upholstery: Artificial leather covered foam, electrically conductive
- Material: Stainless steel
- Power supply: 230 V / Battery: 24 V, 1.2 Ah

Positioning Aids



Head support

Art.-no. 61147

- For positioning of patient in supine position
- Radiolucent
- L x W: 250 x 250mm
- Material: PUR, electrically conductive



Head ring made of silicone gel for adults

Art.-no. 61115

- Height: 45mm | diameter (outer): 200mm

Head ring made of silicone gel for adolescents

Art.-no. 61115_2

- Height: 35mm | diameter (outer): 140mm

Head ring made of silicone gel for infants

Art.-no. 61115_3

- Height: 20mm | diameter (outer): 90mm



Head positioning pad adults made of silicone gel

Art.-no. 61910

- For positioning of patient in prone position
- L x W x H: 280 x 240 x 150mm

Head positioning pad children made of silicone gel

Art.-no. 61911

- For positioning of patient in prone position
- L x W x H: 260 x 230 x 140mm



Head positioning pad

Art.-no. 61905

- For positioning of patient in supine position
- L x W x H: 400 x 500 x 100mm
- Material: Artificial leather covered foam, electrically conductive



Lateral positioning pad

Art.-no. 61915

- For stable positioning of patient during hip, shoulder or kidney surgeries
- Cut-out to relieve pressure from shoulder joint and upper arm
- L x W x H: 700 x 500 x 160 / 25mm
- Material: PU coated foam, electrically conductive



Spinal disc pad

Art.-no. 61900

- For positioning of patient during surgeries at spinal disc
- Radiolucent
- L x W x H: 820 x 450 x 90mm
- Material: Viscoelastic foam, electrically conductive

Positioning Aids



Flat bottom chest roll small made of silicone gel
Art.-no. 61920

- L x W x H: 305 x 100 x 75mm

Flat bottom chest roll medium made of silicone gel
Art.-no. 61921

- L x W x H: 405 x 100 x 75mm

Flat bottom chest roll large made of silicone gel
Art.-no. 61922

- L x W x H: 510 x 100 x 75mm



Bottom roll adults
Art.-no. 61925

- For positioning of patient in supine position
- L x W x H: 500 x 220 x 110mm
- Material: Artificial leather covered foam, electrically conductive

Bottom roll children
Art.-no. 61926

- Same as 61925, but with dimensions
L x W x H: 500 x 150 x 75mm



Foot positioning pad
Art.-no. 61930

- L x W x H: 600 x 220 x 180mm
- Material: Artificial leather covered foam, electrically conductive



Leg positioning tunnel

Art.-no. 61935

- For positioning of patient in lateral position
- Radiolucent
- L x W x H: 600 x 390 x 220mm
- Material: PU coated foam, electrically conductive



Heel pads made of silicone gel (pair)

Art.-no. 61940

- For positioning of patient in supine position
- L x W x H: 190 x 100 x 70mm

Accessories for Anteris™ Ultra 700 Series



Foot switch for Anteris Ultra 700 Series

Art.-no. 704060

- Additional control unit to adjust the main functions height, lateral and Trendelenburg movement
- Tethered



Infrared remote control for Anteris Ultra 700 Series

Art.-no. 704070

- Wireless data transfer
- Control panel for all electric functions of operating table
- Incl. charging station



Longitudinal sliding for Anteris Ultra 700 Series

Art.-no. 71640

- Electrical longitudinal sliding of table top up to 380mm
- Changes lowest height of operating table from 580mm to 600mm
- Cannot be retrofitted!



Motorised drive

Art.-no. 70645

- Motorised fifth castor for manoeuvring and positioning of operating table
- Cannot be retrofitted!



Hand grips for motorised drive

Art.-no. 71570

- For manoeuvring and positioning of operating table by means of electronic drive at operating table's head end
- With integrated speed control
- Material: Stainless steel
- Suitable for: 70645



Back-up system for Anteris Ultra 700 Series

Art.-no. 70655

- System ensures the operability of operating table's basic functions (height, lateral, Trendelenburg and brakes) in case of electronic defect by means of 24V battery
- Cannot be retrofitted!

Accessories for Anteris™ Ultra 700 Series



Radiographic cassette tunnel for Au770
Art.-no. 71642_2

- To insert radiographic cassette carrier

Radiographic cassette tunnel for Au790
Art.-no. 71642_3

- To insert radiographic cassette carrier



Coupling point for extension device 63000_2
Art.-no. 63005_1

- For attachment of extension device, for operating tables of series 700 without sliding

Coupling point for extension device 63000_2
Art.-no. 63005_2

- Same as 63005_1, but for operating tables of series 700 with sliding



Head rest for Anteris Ultra 700 Series
Art.-no. 71144

- Easy adjustment by means of gas spring from -45° to $+45^{\circ}$
- Docking and undocking at back rest
- Radiolucent
- L x W: 266 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable



Head rest for Anteris Ultra 700 Series
Art.-no. 71146

- Easy adjustment by means of gas spring from -45° to $+45^{\circ}$
- Double joint enables additional tilting
- Docking and undocking at back rest
- Radiolucent
- L x W: 266 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable



Head rest for Anteris Ultra 700 Series
Art.-no. 71147

- Easy adjustment by means of gas spring from -45° to $+45^{\circ}$
- Double joint enables additional tilting
- Docking and undocking at back rest
- Radiolucent
- L x W: 266 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable



Head-end rail (piece)
Art.-no. 71148

- Additional rail at head rest
- Length: 300mm
- Material: Stainless steel
- Suitable for: 71144, 71146

Accessories for Anteris™ Ultra 700 Series



Side rails at head rest (pair)

Art.-no. 71149

- Additional side rails at both sides of head rest
- Length: 160mm
- Material: Stainless steel
- Suitable for: 71144, 71146



Upper back rest for Anteris Ultra 700 Series

Art.-no. 71250_2

- Docking and undocking at lower back rest
- Radiolucent
- L x W: 268 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)



Tripartite back rest for Anteris Ultra 700 Series

Art.-no. 71390

- Incl. head rest
- Easy access to shoulder area by means of three detachable segments
- Docking and undocking at lower back rest
- Radiolucent
- L x W: 320 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable



**Pelvis section extension for Anteris Ultra 700 Series
Art.-no. 71125_2**

- Extension of pelvis section by 400mm
- With cut-out
- Radiolucent
- L x W: 400 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)



**Leg rest for Anteris Ultra 700 Series
Art.-no. 71625**

- Electrically adjustable
- Docking and undocking at pelvis section
- Radiolucent
- L x W: 660 x 550mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)



**Divided leg rest for Anteris Ultra 700 Series (pair)
Art.-no. 71615_2**

- Electrically adjustable
- Swivelling up to 180°
- Docking and undocking at pelvis section
- Radiolucent
- L x W: each 660 x 266mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)

Accessories for Anteris™ Ultra 700 Series



Leg rest with access extension for Anteris Ultra 700 Series (pair)

Art.-no. 71620

- Electrically adjustable
- Swivels out for direct patient access
- Docking and undocking at pelvis section
- Radiolucent
- L x W: each 740 x 266mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)



Double divided leg rest for Anteris Ultra 700 Series (pair)

Art.-no. 71660

- Electrically adjustable
- Swivelling up to 180° and bendable
- Docking and undocking at pelvis section
- Radiolucent
- L x W: each 662 x 266mm
- Viscoelastic SF upholstery, H = 80mm, electrically conductive, detachable
- Integrated side rails (EU)



Carbon positioning rest for Anteris Ultra 700 Series


Art.-no. 71630

- Incl. upholstery
- For intraoperative x-ray, optimal radiolucency (360°)
- Docking and undocking at pelvis section or lower back rest
- L x W: 750 x 540mm
- Viscoelastic SF upholstery, H = 70mm, electrically conductive, detachable



**Transfer leg rest for Anteris Ultra 700 Series
Art.-no. 71635**

- For short-time positioning of legs
- Docking and undocking at pelvis section
- Radiolucent
- L x W: 900 x 500mm



**Extension segments (pair)
Art.-no. 61355SF**

- For extension of the operating table's lying surface
- Integrated side rails (EU)
- L x W: each 310 x 110mm
- Viscoelastic SF upholstery, H = 70mm, electrically conductive, detachable
- Suitable for all operating tables with viscoelastic foam upholstery



**Radiographic cassette carrier for Anteris Ultra 700 Series
Art.-no. 71643**

- To insert radiographic cassettes under table top of operating table
- With telescopic rod
- L x W x H: 470 x 385 x 16mm
- Suitable for radiographic cassettes with dimensions 14" x 17" (US)
- Material: Stainless steel
- Suitable for: 71642_2, 71642_3

Accessories for Anteris™ Ultra 600 Series




Upper back rest for models Au680 (L= 360mm)
Art.-no. 61250SF

- Docking and undocking at lower back rest
- Integrated side rails (EU)
- L x W: 360 x 540mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61250)



Upper back rest for models Au680 (L= 250mm)
Art.-no. 61255SF

- Docking and undocking at lower back rest
- Integrated side rails (EU)
- L x W: 250 x 540mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable



Tripartite back rest for models Au680
Art.-no. 61390SF

- Incl. head rest
- Easy access to shoulder area by means of three detachable segments
- Docking and undocking at lower back rest
- L x W: each segment 315 x 170mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61390)



**Stabilization cylinder for Anteris Ultra 600 Series
Art.-no. 50650**

- Extendable
- For optimal stability of operating table
- Cannot be retrofitted!



**Manual back-up system for Anteris Ultra 600 Series
Art.-no. 50655**

- Secures the functionality of operating table in case of an electronic defect or breakdown of power supply
- Override of electronic functions by means of manual foot pedal
- Cannot be retrofitted!

Accessories for Anteris™ Ultra 600 Series



Foot switch for Anteris Ultra 600 Series

Art.-no. 604060

- Additional control unit to adjust the main functions height, lateral and Trendelenburg movement
- Tethered
- Cannot be retrofitted!



Infrared remote control for Anteris Ultra 600 Series

Art.-no. 604070

- Wireless data transfer
- Control panel for all electric functions of operating table
- Incl. charging station



Connecting bar for Anteris Ultra 600 Series

Art.-no. 61670

- Exchange of stainless steel connecting bar at back rest with a radiolucent connecting bar made of special plastic
- Improves radiolucency of back rest



Head rest large for Anteris Ultra 600 Series

Art.-no. 61144SF

- Easy adjustment by means of gas spring from -45° to $+45^{\circ}$
- Docking and undocking at back rest
- L x W: 285 x 495mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61144)



Head rest large with double joint for Anteris Ultra 600 Series

Art.-no. 61146SF

- Easy adjustment by means of gas spring from -45° to $+45^{\circ}$
- Double joint enables additional tilting
- Docking and undocking at back rest
- L x W: 285 x 495mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61146)

Accessories for Anteris™ Ultra 600 Series



**Pelvis section extension for Anteris Ultra 600 Series
(L= 430mm)**

Art.-no. 61125SF

- Extension of pelvis section by 430mm
- With cut-out
- Integrated side rails (EU)
- L x W: 430 x 540mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61125)



**Pelvis section extension for Anteris Ultra 600 Series
(L= 250mm)**

Art.-no. 61130SF

- Extension of pelvis section by 250mm
- With cut-out
- Integrated side rails (EU)
- L x W: 250 x 540mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable



Leg rest for Au680

Art.-no. 61625SF

- Adjustment by means of gas spring from -90° up to $+30^{\circ}$
- Docking and undocking at pelvis section
- Integrated side rails (EU)
- L x W: 670 x 540mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61625)



Divided leg rest for Au680 (pair)

Art.-no. 61615SF

- Adjustment by means of gas spring from -90° up to $+30^{\circ}$, swivelling up to 180°
- Docking and undocking at pelvis section
- Integrated side rails (EU)
- L x W: each 730 x 250mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61615)



Double divided leg rest for Au680 (pair)

Art.-no. 61660SF

- Adjustment by means of gas spring from $+30^{\circ}$ up to -90° , swivelling up to 180° , bendable
- Docking and undocking at pelvis section
- Integrated side rails (EU)
- L x W: each 720 x 250mm
- Upholstery: Viscoelastic SF upholstery, H = 80mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61660)

Accessories for Anteris™ Ultra 600 Series



Carbon positioning rest for Anteris Ultra 600 Series
Art.-no. 61690

- Incl. upholstery
- For intraoperative x-ray, optimal radiolucency (360°)
- Docking and undocking at pelvis section or lower back rest
- L x W: 750 x 540mm
- Upholstery: SF, H = 45mm, radiolucent, electrically conductive



Paediatric surgery rest for Anteris Ultra 600 Series
Art.-no. 61633

- Docking and undocking at pelvis section
- Integrated side rails (EU)
- L x W: 700 x 430mm
- Upholstery: Viscoelastic foam, H = 60mm, electrically conductive



Transfer leg rest for Anteris Ultra 600 Series
Art.-no. 61635

- For short-time positioning of legs
- Docking and undocking at pelvis section
- L x W: 900 x 500mm



Extension segments all models (pair)

Art.-no. 61355SF

- Segments for extension of operating table's lying surface
- Integrated side rails (EU)
- L x W: each segment 310 x 110mm
- Upholstery: Viscoelastic SF upholstery, H = 70mm, radiolucent, electrically conductive, detachable
- Also available with PUR upholstery, H = 60mm (Art.-no. 61355)



Rinsing set for Anteris Ultra 600 Series

Art.-no. 61120

- Incl. bucket, basin, strainer, tube, undercarriage and fixation
- Docking and undocking at pelvis section
- L x W x H (Basin): 325 x 265 x 65mm
- Volume (Bucket): 15 litre
- Material: Stainless steel



Radiographic cassette carrier for Anteris Ultra 600 Series

Art.-no. 61640

- To insert radiographic cassettes under table top of operating table
- With telescopic rod
- L x W x H: 470 x 385 x 16mm
- Suitable for radiographic cassettes with dimensions 14" x 17" (US)
- Material: Stainless steel

General Accessories

61511	Rotatable clamp	4	61351	Shoulder support (pair)	9
61512	Rotatable clamp with quick lock and release lever	4	61352	Shoulder support incl. concave upholstery (pair)	9
61551	Clamp for square rod (25 x 10mm)	4	61353	Shoulder support incl. straight upholstery (pair)	10
61552	Clamp for square rod (20 x 20mm)	5	61362	Double side support	10
61553	Clamp for round rod (18mm)	5	61367	Joint arm for body supports	10
61554	Clamp with quick lock and release lever for round rod (18mm)	5	61361	Holder for body supports	11
61750	Anaesthesia screen	6	61361_1	Body support for pubis, sacrum and sternum	11
61755	Arm restraint cuff	6	61361_2	Body support for side and shoulder	11
61751	Flexible anaesthesia screen	6	61361_3	Body support for back and buttocks	12
61752	Tube support	7	61361_4	Supporting roll	12
61800	Infusion holder	7	91611	Instrument disposal round form	12
61805	Rotatable infusion holder	7	91612	Instrument disposal squared form	12
61410	Mayo table	8	91613	Bowl with drainage and tube	13
61560	Side rail extension	8	91608	Bucket	13
61571	Sliding handles (pair)	8	61575	Hand switch support	13
61341	Side support (pair)	9			

Head Positioning

61100	Universal adapter for Anteris Ultra 600 Series	14
71100	Universal adapter for Anteris Ultra 700 Series	14
61114_1	Horseshoe shaped head rest with cut-out	14
61114_3	Horseshoe shaped head rest two-parted	14
61395	Head support	15
61117	Joint arm	15
61117_2	Head support for joint arm	15
61117_3	Horseshoe shaped head rest with cut-out for joint arm	16
61117_5	Horseshoe shaped head rest two-parted for joint arm	16
61113	Wrist support ring for head positioning accessories	16

Arm Positioning

61211	Arm rest	18
61227	Beach-Chair- shoulder OP arm positioner	18
61225	Arm rest with joint arm and one-hand control	18
61222	Arm rest with one-hand control and comfort upholstery (Length 600mm)	19
61223	Arm rest with one-hand control and comfort upholstery (Length 450mm)	19
61220SF	Arm rest standard	19
61092	Side support / arm rest	20
71345	Arm protection shield made of stainless steel	20
61400	Hand cuffs (pair)	20

Leg Positioning

61700	Leg rest according to Goepel (piece)	22
81705	Leg rest according to Goepel Comfort (pair)	22
61725	Stirrups with one-hand-adjustment, gasspring assisted (pair)	22
61727	Secure-Release™ Clamp, clip-on, EU	23
61728	Select-Blade Clamp, suspendable, EU	23
61710	Leg rest for adiposity patients (pair)	23
61780	Elbow support (pair)	24
61760	Lower leg fixation	24
61770	Thigh fixation	24
61632	Rectal positioning device for Anteris Ultra 600 Series	25
71632	Rectal positioning device for Anteris Ultra 700 Series	25
61150	Knee supports for Anteris Ultra 600 Series (piece)	25
61315	Rectoscopy roll	26
61630	Knee support	26
61631	Foot plates (pair)	26
61405	Body and leg strap	27
61406	Leg strap (pair)	27
61407	Leg strap (pair)	27

Neurosurgery

661065	Standard holder	28
661067	Swivel adapter	28
661069	Skull clamp	28
661073	Skull pins adults	29
661074	Skull pins children	29
661071	Cross bar adapter	29
661076	Horseshoe headrest for adults	29
661078	Horseshoe headrest for children	29

Orthopaedics and Traumatology

63000_1	Extension device for Anteris Ultra 600 Series	30
63000_2	Extension device for Anteris Ultra 700 Series	30
63005_1	Coupling point for extension device 63000_2 for Anteris Ultra 700 Series operating tables without sliding	31
63005_2	Coupling point for extension device 63000_2 for Anteris Ultra 700 Series operating tables of with sliding	31
61010	Extension shoe for adults (piece)	31
61020	Extension shoe for children	31
63025_1	Countertraction post for femur for extension device 63000_1	31
63025_2	Countertraction post for femur for extension device 63000_2	31
63045	Countertraction post for tibia and fibula	32
63030	Pelvis plaster cast unit	32
63060	Horseshoe holder "Kirschner"	32
63026	Leg support	33
63070_1	Transfer leg rest (pair)	33
63070_2	Transfer leg rest	33
63090	Trolley for extension device	33
64000	Radiolucent extension device for Anteris Ultra 600 Series	34
64010	Extension shoe	34
62050	Arm extension device	35
62055	Weinberger hand cuff	35
61226	Upper arm support	35
61231	Arm and hand surgical table levitating	36
61232	Arm and hand surgical table with supporting leg	36
81230	Arm and hand surgical table made of carbon	36
61330	Meniscotomy roll	37
61328	Arthroscopic leg holder system	37
61327	Knee positioning device	37

Positioning Aids

61147	Head support	38
61115	Head ring made of silicone gel for adults	38
61115_2	Head ring made of silicone gel for adolescents	38
61115_3	Head ring made of silicone gel for infants	38
61910	Head positioning pad adults made of silicone gel	38
61911	Head positioning pad children made of silicone gel	38
61905	Head positioning pad	39
61915	Lateral positioning pad	39
61900	Spinal disc pad	39
61920	Flat bottom chest roll small made of silicone gel	40
61921	Flat bottom chest roll medium made of silicone gel	40
61922	Flat bottom chest roll large made of silicone gel	40
61925	Bottom roll adults	40
61926	Bottom roll children	40
61930	Foot positioning pad	40
61935	Leg positioning tunnel	41
61940	Heel pads made of silicone gel (pair)	41

Accessories for Anteris™ Ultra 700 Series

704060	Foot switch for Anteris Ultra 700 Series	42	71390	Tripartite back rest for Anteris Ultra 700 Series	46
704070	Infrared remote control for Anteris Ultra 700 Series	42	71125_2	Pelvis section extension for Anteris Ultra 700 Series	47
71640	Longitudinal sliding for Anteris Ultra 700 Series	42	71625	Leg rest for Anteris Ultra 700 Series	47
70645	Motorised drive	43	71615_2	Divided leg rest for Anteris Ultra 700 Series (pair)	47
71570	Hand grips for motorised drive	43	71620	Leg rest with access extension for Anteris Ultra 700 Series (pair)	48
70655	Back-up system for Anteris Ultra 700 Series	43	71660	Double divided leg rest for Anteris Ultra 700 Series (pair)	48
71642_2	Radiographic cassette tunnel for Au770	44	71630	Carbon positioning rest for Anteris Ultra 700 Series	48
71642_3	Radiographic cassette tunnel for Au790	44	71635	Transfer leg rest for Anteris Ultra 700 Series	49
63005_1	Coupling point for extension device 63000_2	44	61355SF	Extension segments (pair)	49
63005_2	Coupling point for extension device 63000_2	44	71643	Radiographic cassette carrier for Anteris Ultra 700 Series	49
71146	Head rest for Anteris Ultra 700 Series	45			
71147	Head rest for Anteris Ultra 700 Series	45			
71148	Head-end rail (piece)	45			
71149	Side rails at head rest (pair)	46			
71250_2	Upper back rest for Anteris Ultra 700 Series	46			

Accessories for Anteris™ Ultra 600 Series

61250SF	Upper back rest for models Au680 (L= 360mm)	50
61255SF	Upper back rest for models Au680 (L= 250mm)	50
61390SF	Tripartite back rest for models Au680	50
50650	Stabilization cylinder for Anteris Ultra 600 Series	51
50655	Manual back-up system for Anteris Ultra 600 Series	51
604060	Foot switch for Anteris Ultra 600 Series	52
604070	Infrared remote control for Anteris Ultra 600 Series	52
61670	Connecting bar for Anteris Ultra 600 Series	52
61144SF	Head rest large for Anteris Ultra 600 Series	53
61146SF	Head rest large with double joint for Anteris Ultra 600 Series	53
61125SF	Pelvis section extension for Anteris Ultra 600 Series (L= 430mm)	54
61130SF	Pelvis section extension for Anteris Ultra 600 Series (L= 250mm)	54
61625SF	Leg rest for Ar680	55
61615SF	Divided leg rest for Ar680 (pair)	55
61660SF	Double divided leg rest for Ar680 (pair)	55
61690	Carbon positioning rest for Anteris Ultra 600 Series	56
61633	Paediatric surgery rest for Anteris Ultra 600 Series	56
61635	Transfer leg rest for Anteris Ultra 600 Series	56
61355SF	Extension segments all models (pair)	57
61120	Rinsing set for Anteris Ultra 600 Series	57
61640	Radiographic cassette carrier for Anteris Ultra 600 Series	57



United Kingdom



EMEA



Asia Pacific



Latin & South America



South Asia



USA & Canada



Contact the team
+44 (0)113 2777 393
enquiries@brandon-medical.com
www.brandon-medical.com



In line with the Brandon policy of progressive improvement, the right is reserved to alter specifications and/or prices without notice. Brandon Medical, Quasar, HD-LED and Symposia are registered trademarks of Brandon Medical Co Ltd. Anteris Accessories A1-20525 20MAY2025