

# P2 Series

# **General Product Overview**

Brandon Medical Ceiling Mounted Retractable Pendants are supplied at the customer specified height. The 300x300mm square column is powder coat finished to RAL9010 and supplied with a ceiling cover. Every pendant is bespoke, allowing for a wide range of 13A switched sockets, data sockets, equipotential earth sockets as well as up to 9 Medical Gas Terminals.

Height can be adjusted by up to 300mm using the electrically powered motor.

#### **Main Features**

- Modular and configurable design
- P2 series is available in different lengths to suit specific requirements
- Vertical mounted gas terminals. All terminal units used for medical gas are: Penlon East / Medaes Gem 10
- Low maintenance
- Manufactured in the UK


### Standards

- HTM 02-01: Health Technical Memorandum Medical Gas Pipeline Systems
- BS EN 7396: Medical Gas Pipelines for Compressed Medical Gases and Vacuum
- BS EN 9170:
- Terminal Units for use with Compressed Medical Gases and Vacuum
- BS EN 15908:

Non-interchangeable screw Threaded (NIST) Low-pressure Connectors for Medical Gases









1/ /	
Item 1	Retractable Pendant Body
Item 2	Pendant Outer Body
Item 3	Ceiling Cover
Item 4	Adjustable Ceiling Spacer
Item 5	Adjustable Drop Rods

Description	Weight
First Fix Spacer	31kg
Pendant	70kg
Total	101kg





Soffit Plate





ELMFIELD ROAD | MORLEY | LEEDS | LS27 OEL | WEST YORKSHIRE | ENGLAND. T: +44 113 277 7393 / 0845 1243 666 F: +44 113 272 8844 / 0845 1243 667 E: enquires@brandon-medical.com W: www.brandon-medical.com Brandon Medical Company Limited. Registered in England & Wales No. 2827189. VAT No. GB734 591421. A subsidiary of Brandon Group Limited.