


Quasar® eLite Mobile High Performance Operating Light



Near Perfect Colour Rendition with Red Balance Control

For the best visualisation of tissue during surgery.

Fat Beam Technology

The largest light beams of any current product with high intensity across the full beam width.

Unique Comfort Halo

A soft halo of light to reduce eye strain from high contrasts of light intensity.



Designed to Minimise Infections

Easy to clean, sealed light head, remote controls and anti-microbial coatings.

Premium Quality Movement

Low weight light head and perfectly balanced arm system.



Light Quality Characteristics

Maximum light intensity 160,000 Lux

Beam size characteristics
(measured at 1m)

Light field diameter (d10) 200-420mm

Light field diameter (d50) 120-250mm

Beam uniformity (d50/d10
ratio) 0.6



Electrical & Mechanical Characteristics

Nominal power consumption 48W

Radiant energy mWm².Lux <3.31

IP rating of light engine 54

LED life > 60,000 hours

Part Number

Qe60m Quasar® eLite mobile

Qe60ms Quasar® eLite mobile Schuko style plug

Qe60mu Quasar® eLite mobile US style plug

Qe60mel Quasar® eLite mobile with 3 hr emergency power

Brandon Medical Co Ltd, Elmfield Road, Morley, Leeds, LS27 0EL

T: +44 (0)845 1243 666 F: +44 (0)845 1243 667

T: +44 (0)113 2777 393 F: +44 (0)113 2728 844

E: enquiries@brandon-medical.com

www.brandon-medical.com